

4th June, 2018

- Incubator has been confirmed sterile.
- Reviving HEK293 cells from frozen stock.
- ↳ ~~All~~ All three micropipettes, pipette fillers, 2 boxes of 1mL tips, 3 boxes 200µL tips (both filters), and about five or six Falcon 75mL tubes have been left permanently in the hood [Falcon tubes borrowed from Arankshet].
- ↳ ~~About~~ One vial of frozen HEK293 stock from SG's lab, obtained through Rutika - Added to 9mL of thawed DMEM [aliquot given by Rutika], spun at 2000rpm for 10 minutes to pellet out cells.
- ↳ Cell pellet, ^{has been} resuspended in 1mL of DMEM, and 500µL each has been ~~to~~ added drop-wise to 2 60mm dishes with 5mL DMEM.
- ↳ ⇒ 2 60mm dishes seeded with HEK293 cells of passage #6.

