

6th June, 2018:

→ Passaging of P6 cells-

- ↳ Cells are at 100% confluency
- ↳ One 60mm plate trypsinized first - Trypsin wash followed by treatment - both with 500µL 0.25% trypsin - EDTA [aliquot from Rutika].
- ↳ After 1 min incubation, trypsin was inhibited with an equal amount of DMEM-FBS - we ^{have} used about 5mL of the mix aliquot as we only have 15mL Falcon tubes, but ^{5mL} Eppendorf tubes would be ideal for these volumes (<1mL).
- ↳ Cells are ~~the pellet~~ ~~pelleted~~ spun down in about 6mL of DMEM:FBS (pooled from both 60mm dishes).
- ↳ This pooled set has split between 2 100mm dishes and one frozen stock.
- ↳ Proportion of components are roughly 1:2 from 60mm to 100mm plates - about 12-15mL of media and ~1mL of cell suspension.



- ↳ In this batch, pen-strep [aliquot from Rutika] was added directly to the dishes as we forgot to add to the DMEM: FBS aliquot - 200µL in each dish, and in the remaining media aliquot.
 - ↳ Proportions for 1mL stock - 700µL media + 200µL cells + 100µL DMSO.
 - ↳ Finally: 2 100mm dishes of ~~100mm~~ HEK293 - P7
 & 1 cryovial of 1mL stock of HEK293-P6.
- No passaging is expected tomorrow.