

30th May, 2018

- Fogging of the culture room
 - ↳ Fogging machine borrowed from ML Lab through Aishwarya
 - ↳ 20ml of fogging solution gotten through Kalpesh
 - ↳ 3% solution in MilliQ ~~water~~ (from 2nd floor fly lab), of which a little more than 1L was used. 10ml of concentrate remaining with us [10ml in ~990ml of water].
 - ↳ Started at 10:30 PM, after covering up the microscope with as much tape and plastic covers as possible.
 - ↳ Incubator trays cleaned with soap and kept outside.
 - ↳ Both the laminar hood and the incubator have been open, as have been the one set of shelves.
 - ↳ Run has been terminated at 11:30 PM.
- A few instruction manuals and spare serological pipettes have been found in the room, and stored in the top locked shelf (key #0063)
- Key for the room is yet to be issued to the lab
- The second laminar hood, with button controls for the glass separation, has been left closed as it isn't going to be used.